



Product List 2012



Content

Luminex

Luminex System	1
Multimetrix Assays	1
AtheNA Multi-Lyte Assays	1
Signature Assays	2
Accessories for Luminex System	2

Molecular Biology

Illumigene System	4
Asuragen PCR Test	4
Euroclone PCR	5
Prodesse Realtime PCR	7
PCR Fast Tests	8

Rapid Tests

Virology	10
Bacteriology	10
Parasitology	11
Clinical Chemistry	12

ELISA - Antibody Detection

Virology	13
Bacteriology	16
Parasitology	18
Mycology	19
Immunology	19

ELISA - Antigen Detection

Virology	21
Bacteriology	21
Parasitology	21
Mycology	21
Clinical Chemistry	22

Immunoblot

Virology	22
Immunology	22

IFA - Antibody Detection

Virology	23
Bacteriology	24
Parasitology	25
Immunology	25

IFA - Antigen Detection

Virology	26
Bacteriology	26
Parasitology	26
Mycology	26

Agglutination Tests

Virology	27
Bacteriology	27
Mycology	27

Collection and Processing Systems for Stool

Closed Collection & Concentration	28
Collection Systems with Integrated	28
Filtering Centrifuge Tube	28
Reagents	29

Other

Mycology	29
----------	----

Reagents

Reagents	30
----------	----

Luminex System

Luminex	LX200 IS	Luminex 200 IS - System	Analyzer, XY-Platform, Sheath Delivery System, PC, Flatscreen 17", Inkjet Printer, Software v 2.3, Barcodescanner, 1 Set Calibrators / Controls, 20L Sheath Fluid.
---------	----------	--------------------------------	--

Multimetrix Assays

Borrelia				
Progen Biotechnik	MM2001 201004	Borrelia IgM	Multiplex test for the detection of IgM against a Borrelia lysate and 8 recombinant Borrelia afzelii/garinii antigens (OspA, OspC (2x), p100, p18, p39, p41-I, VlsE-C6)	96 tests
Progen Biotechnik	MM2002 201003	Borrelia IgG	Multiplex test for the detection of IgG against a Borrelia lysate and 8 recombinant Borrelia afzelii/garinii antigens (OspA, OspC, p100, p18, p39, p41-I, p58, VlsE-C6)	96 tests
Progen Biotechnik	MM2014	Borrelia IgM Screen	Multiplex test for the detection of IgM against a lysate of Borrelia afzelii and recombinant Osp C of Borrelia afzelii and garinii	96 tests
Progen Biotechnik	MM2015	Borrelia IgG Screen	Multiplex test for the detection of IgG against 2 lysate antigens (Borrelia afzelii and garinii)	96 tests
HPV				
Progen Biotechnik	MMLK0124	HPV Genotyping	Multiplex test for research in epidemiology on 24 Human Papillomaviruses (HPV) in PCR amplified samples, according to Schmitt et al. 2006	96 tests

AtheNA Multi-Lyte Assays

Autoimmune Diagnostics				
Zeus Scientific	A91001M 201034	RF Plus	Test system for the detection of Rheumafactor antibodies	96 tests
Zeus Scientific	A21101 201023	ANA II Plus	Multiplex test for the detection of 10 different autoimmune antibodies against SSA, SSB, Sm, RNP, Scl-70, Jo-1, Centromere B, Histones, dsDNA and HEp2NA	96 tests
Zeus Scientific	A22101 201024	ANA III Plus	Multiplex test for the detection of 10 different autoimmune antibodies against SSA-52, SSA-60, SSB, Sm, RNP, Scl-70, Jo-1, Centromere B, Ribosomal P, dsDNA	96 tests
Zeus Scientific	A96101 201025	Autoimmune Vasculitis Plus	Multiplex test for the detection of autoimmune antibodies against MPO, PR3 and GBM	96 tests
Zeus Scientific	A51101 201017	Tg/TPO Plus	Multiplex test for the detection of antibodies against Tg and TPO	96 tests
Zeus Scientific	A61101G 201015	Celiac IgG Plus	Multiplex test for the detection of IgG against Gliadin and human tissue transglutaminase	96 tests
Zeus Scientific	A61101A 201014	Celiac IgA Plus	Multiplex test for the detection of total IgA and IgA against Gliadin and human tissue transglutaminase	96 tests

AtheNA Multi-Lyte Assays

EBV				
Zeus Scientific	A92101G 201032	EBV IgG Plus	Multiplex test for the detection of IgG against VCA, EBNA and EA	96 tests
Zeus Scientific	A92201M 201031	EBV IgM Plus	Multiplex test for the detection of IgM against VCA	96 tests
HSV				
Zeus Scientific	A97101G	HSV 1&2 IgG Plus	Multiplex test for the detection of IgG against HSV-1 and HSV-2 (recombinant gG-1 and gG-2)	96 tests
Measles, Mumps, Varizella				
Zeus Scientific	A93111G	MMV IgG Plus	Multiplex test for the detection of IgG against Measles virus, Mumps virus and VZV	96 tests
Zeus Scientific	A93101G	MMRV IgG Plus	Multiplex test for the detection of IgG against Measles, Mumps, Rubella virus and VZV	96 tests
Toxoplasma, Rubella, CMV, HSV				
Zeus Scientific	A81101G 201018	ToRCH IgG Plus	Multiplex test for the detection of IgG against Toxoplasma, Rubella, CMV, HSV-1 and HSV-2	96 tests
Zeus Scientific	A81101M 201019	ToRCH IgM Plus	Multiplex test for the detection of IgM against Toxoplasma, Rubella, CMV, HSV-1 and HSV-2	96 tests

Signature Assays

KRAS/BRAF				
Asuragen	49418	KRAS/BRAF Mutations	Multiplex assay for detection of 12 mutations in codons 12 and 13 of KRAS gene and BRAF mutation V600E	48 rxn
Leukemia				
Asuragen	76004	LTx	Leukemia Translocation Panel. Multiplex assay for detection of 12 fusion transcripts from 9 translocations.	24 rxn
Asuragen	46072	NPM1 Mutations	Multiplex assay for detection of transcripts A, B and D of NPM1 gene (RUO)	24 rxn

Accessories for Luminex System

Accessories Luminex				
Luminex	MM1101 201081	Sheath Fluid		20 L
Luminex	MM1112 201080	Sheath Fluid concentrate (20x)		1 L
Luminex	MM1102 201082	Calibrator Cal1		5 ml

Accessories for Luminex System

Accessories Luminex				
Luminex	MM1103 201083	Calibrator Cal2		5 ml
Luminex	MM1104 201084	Control Con1		5 ml
Luminex	MM1105 201085	Control Con2		5 ml
	MM1111 201091	Microtiterplate C8 Polysorb		1 pce
Zeus Scientific	101D 368001	AtheNA Multi-Lyte 96 Well Filterplate		1 pce
Zeus Scientific	AX9001 201094	AtheNA Multi-Lyte Daily Maint. Reagents	provides all of the materials necessary to perform 24 daily startup and 24 daily shutdown procedures on the Luminex system	
Zeus Scientific	005M 201098	Sample Diluent	for EBV samples with INV:NSC	30 ml
	MM1006 201092	Multiscreen HTS Vacuum Manifold		1 pce
	MM1011 201095	BioTek Microplate Strip Washer	For automated washing of filterplates (Model ELx50/8FMW)	1 pce
	MM1003 201093	Ultrasonics Bath Sonorex RK 31		1 pce

Molecular Biology

illumigene System

Meridian Bioscience		illumigene	<i>illumigene is a simplified new technology platform, based upon Loop Mediated Isothermal Amplification (LAMP) that makes molecular biology testing possible for any laboratory. The very compact illumipro-10 incubator/reader can run up to 10 specimens and results are available in less than one hour.</i>	
Meridian Bioscience	280050 308520	illumigene C. difficile	DNA amplification assay for the detection of toxigenic Clostridium difficile	50 tests
Meridian Bioscience	279920 308524	illumigene C. difficile External Control	Positive and Negative Control Reagents for use with the illumigene C. difficile kit	18 tests
Meridian Bioscience	280350	illumigene GBS	DNA Amplification Assay for the Detection of Group B Streptococcus in Vaginal/Rectal Antepartum specimens.	50 tests
Meridian Bioscience	279900	illumigene GBS External Control	Positive and Negative Control Reagents for use with the illumigene GBS kit	18 tests
Meridian Bioscience	610172 308501	illumipro-10	Automated isothermal amplification and detection system	1 pce

Asuragen PCR Tests

Fragile X				
Asuragen	76008 380103	AmplideX FMR1	Two PCR technologies for molecular characterization of CGG trinucleotide repeat in 5' primer region of FMR1 gene	100 rxn
Asuragen	49442	AmplideX FMR1 mPCR	Innovative PCR only approach for the detection of methylation status in the FMR1 gene and elimination of the need for Southern Blot. Amplifies pre-mutation and full mutation alleles up to 1300 or more CGG repeats.	24 rxn
KRAS/BRAF				
Asuragen	49418	KRAS/BRAF Mutations	Multiplex assay for detection of 12 mutations in codons 12 and 13 of KRAS gene and BRAF mutation V600E	48 rxn
Leukemia				
Asuragen	76004	LTx	Leukemia Translocation Panel. Multiplex assay for detection of 12 fusion transcripts from 9 translocations.	24 rxn
Asuragen	76005	BCR/ABL1 Quant	Quantitative detection of BCR/ABL1 fusion transcripts (b2a2, b3a2, and e1a2) using Taqman technology.	48 rxn
Asuragen	49462	ARQ IS Calibrator Panel (e13a2)	Secondary Reference Material for the Standardization of BCR-ABL1 Quantitative Measurement on the International Scale (IS). 4 levels of e13a2; negative control.	1 set
Asuragen	49463	ARQ IS Calibrator Panel (e14a2)	Secondary Reference Material for the Standardization of BCR-ABL1 Quantitative Measurement on the International Scale (IS). 4 levels of e14a2; negative control.	1 set
Asuragen	46072	NPM1 Mutations	Multiplex assay for detection of transcripts A, B and D of NPM1 gene (RUO)	24 rxn

Euroclone PCR Tests

Cardiovascular diseases				
Euroclone	EER001032	Factor II G20210A Genotyping	Duplica Real Time allelic discrimination kit. Provided with two ready-to-use reaction mixes (Amplification Mix and Oligo Mix).	32 tests
Euroclone	EER002032	Factor V G1691A Genotyping	Duplica Real Time allelic discrimination kit. Provided with two ready-to-use reaction mixes (Amplification Mix and Oligo Mix).	32 tests
Euroclone	EER010032	Factor V 1299 Genotyping	Duplica Real Time allelic discrimination kit. Provided with two ready-to-use reaction mixes (Amplification Mix and Oligo Mix).	32 tests
Euroclone	EER013032	Factor V Y1702C Genotyping	Duplica Real Time allelic discrimination kit. Provided with two ready-to-use reaction mixes (Amplification Mix and Oligo Mix).	32 tests
Euroclone	EER024032	Factor V R306T (Cambridge)	Duplica Real Time allelic discrimination kit. Provided with two ready-to-use reaction mixes (Amplification Mix and Oligo Mix).	32 tests
Euroclone	EER022032	Factor XIII Genotyping	Duplica Real Time allelic discrimination kit. Provided with two ready-to-use reaction mixes (Amplification Mix and Oligo Mix).	32 tests
Euroclone	EER003032	MTHFR C677T Genotyping	Duplica Real Time allelic discrimination kit. Provided with two ready-to-use reaction mixes (Amplification Mix and Oligo Mix).	32 tests
Euroclone	EER007032	MTHFR A1298C Genotyping	Duplica Real Time allelic discrimination kit. Provided with two ready-to-use reaction mixes (Amplification Mix and Oligo Mix).	32 tests
Euroclone	EER008032	MTRR A66G Genotyping	Duplica Real Time allelic discrimination kit. Provided with two ready-to-use reaction mixes (Amplification Mix and Oligo Mix).	32 tests
Euroclone	EER009032	MTR A2756G Genotyping	Duplica Real Time allelic discrimination kit. Provided with two ready-to-use reaction mixes (Amplification Mix and Oligo Mix).	32 tests
Euroclone	EER017032	PAI 1 Genotyping kit	Duplica Real Time allelic discrimination kit. Provided with two ready-to-use reaction mixes (Amplification Mix and Oligo Mix).	32 tests
Euroclone	EES001020	RhyMA test CVD	Multiparametric detection of 5 different polymorphisms in cardiovascular diseases (Factor V R506Q and H1299R, Factor II G20210A, MTHFR C677T and A1298C). Multiplex DNA amplification and hybridization of amplified and biotinylated DNA fragments with strips.	20 tests
Euroclone	EES002020	RhyMA test OMO	Multiparametric detection of 5 different polymorphisms in cardiovascular diseases and cancer (MTHFR C677T and A1298C, MS A2756G, MTRR A66G, DNMT3B C149T). Multiplex DNA amplification and hybridization of amplified and biotinylated DNA fragments with strips.	20 tests
Clostridium difficile				
Euroclone	EBR026032	Clostridium difficile A/B	Duplica Real Time Advanced PCR Kit	32 tests
Genetic diseases				
Euroclone	EER005032	Haemochromatosis C282Y	Duplica Real Time allelic discrimination kit. Provided with two ready-to-use reaction mixes (Amplification Mix and Oligo Mix).	32 tests
Euroclone	EER006032	Haemochromatosis H63D	Duplica Real Time allelic discrimination kit. Provided with two ready-to-use reaction mixes (Amplification Mix and Oligo Mix).	32 tests
Euroclone	EER019032	Haemochromatosis S65C	Duplica Real Time allelic discrimination kit. Provided with two ready-to-use reaction mixes (Amplification Mix and Oligo Mix).	32 tests

Molecular Biology

Euroclone PCR Tests

Genetic diseases				
Euroclone	EES006020	RhyMA test HFE	Reverse hybridisation on strips	20 tests
Euroclone	EER025032	Gilbert's Syndrome	Duplica Real Time PCR Kit	32 tests
HSV				
Euroclone	EES005020	RhyMA Test - Herpes screening	Reverse hybridisation on strips	20 tests
Male fertility disorders				
Euroclone	EDP003024	Chromosome Y Microdeletion	Agarose Gel Detection PCR	40 tests
Euroclone	EES003020	RhyMA Test IVF AZF	Multiparametric detection of the Y Chromosome microdeletions located in regions AZFa, AZFb, AZFc according to the guidelines of EMQN 2001. Multiplex DNA amplification and hybridization of amplified and biotinylated DNA fragments with strips.	20 tests
Pharmacogenetics				
Euroclone	EER021060	HLA-B*5701 Genotyping	Duplica Real Time PCR Kit	60 tests
Euroclone	EES007020	RhyMA Test Warfarin	Reverse hybridisation on strips	20 tests
Respiratory infections				
Euroclone	EBR019032	Mycobacterium tuberculosis	Duplica Real Time PCR Kit	32 tests
Euroclone	EBR007032	Legionella pneumophila	Duplica Real Time PCR Kit	32 tests
Euroclone	EBR008032	Adenovirus	Duplica Real Time PCR Kit	32 tests
Euroclone	EBR010032	Mycobacterium avium complex	Duplica Real Time PCR Kit	32 tests
Euroclone	EBR011032	Legionella ssp	Duplica Real Time PCR Kit	32 tests
Euroclone	EBR012032	Mycoplasma pneumoniae	Duplica Real Time PCR Kit	32 tests
Euroclone	EBR006032	Chlamydia pneumoniae	Duplica Real Time PCR Kit	32 tests
Sexually transmitted diseases				
Euroclone	EBR009032	Trichomonas vaginalis	Duplica Real Time PCR Kit	32 tests
Euroclone	EBR018032	Neisseria gonorrhoeae	Duplica Real Time PCR Kit	32 tests
Euroclone	EBR027032	Mycoplasma genitalium	Duplica Real Time Advanced PCR Kit	32 tests
Euroclone	EBR015032	Mycoplasma hominis	Duplica Real Time PCR Kit	32 tests

Molecular Biology

Euroclone PCR Tests

Sexually transmitted diseases

Euroclone	EBR016032	Ureaplasma urealyticum	Duplica Real Time PCR Kit	32 tests
Euroclone	EBR025032	Chlamydia trachomatis	Duplica Real Time Advanced PCR Kit	32 tests

Prodesse Realtime PCR

Clostridium difficile

Prodesse	E38VK00	ProGastro Cd Detection Kit	Real Time Multiplex PCR Kit for the detection of toxigenic strains of Clostridium difficile (tcdB Gene). Including Supermix (with taq), internal and positive control.	100 tests
Prodesse	E38VK55	ProGastro Cd Buffer Kit	Stool Transport and Recovery Buffer for 200 reactions.	100 ml

hMPV

Prodesse	H27VK00	Pro hMPV+	Real Time Multiplex PCR Kit for the detection of Human Metapneumovirus. Including Supermix (with taq), MMLV, RNase Inhibitor, Internal and Positive Controls.	100 tests
----------	---------	------------------	---	-----------

Influenza Virus

Prodesse	H44VK00	ProFlu+ Detection Kit	Real Time Multiplex PCR Kit for the detection of Influenza A, Influenza B, RSV. Including Supermix (with taq), MMLV, RNase Inhibitor and Internal Control	100 tests
Prodesse	H44VK55	ProFlu+ Control Kit	Positive Controls for ProFlu+	4x300µl
Prodesse	H33CK00	ProFlu-ST	Real Time Multiplex PCR Kit. Subtypes Influenza A – Seasonal H1, Seasonal H3, 2009 H1N1. Including Supermix (with taq), MMLV, RNase Inhibitor, Internal and Positive Control.	100 tests

Mycoplasma pneumoniae

Prodesse	P53SK96	ProPneumo-1	Real Time Multiplex PCR Kit for the detection of C. pneumoniae and M. pneumoniae. Including Supermix, Positive Controls and Internal Control. Does not include enzymes.	96 tests
----------	---------	--------------------	---	----------

Parainfluenza Virus

Prodesse	H81VK00	ProParaflu+	Real Time Multiplex PCR Kit for the detection of Parainfluenza Virus 1,2 and 3. Including Supermix (with taq), MMLV, RNase Inhibitor and Internal Control	100 tests
----------	---------	--------------------	---	-----------

PCR Fast Tests

<i>ifp</i>	PCRFast Test Kits	<p><i>Realtime PCR with hybridization probes and internal controls</i></p> <p><i>All PCRFast Kits contain dried reagents with specific primers, probes (FAM) and internal amplification- or transcription controls (HEX).</i></p> <p><i>2 test formats available:</i></p> <p><i>A. 0.2 ml reaction vials ready-to-use</i> <i>88 reaction vials for samples and extraction- and negative controls.</i> <i>8 reaction vials for positive controls</i> <i>(96 reactions in total)</i> <i>- for block cyclers</i></p> <p><i>B. 0.5 ml vials for aliquots</i> <i>9 vials for reconstitution of 10 reactions each for analysis of samples and extraction- and negative controls. 9 vials for reconstitution of 1 reaction each for positive controls</i> <i>(99 reactions in total)</i> <i>- for all PCR cyclers</i></p>	96 rxn
------------	--------------------------	---	--------

Adenovirus

<i>ifp</i>	<i>PAD0180S</i>	Adenovirus	Qualitative detection of Adenovirus in faeces in realtime	96 rxn
------------	-----------------	-------------------	---	--------

Astrovirus

<i>ifp</i>	<i>PAST0190S</i>	Astrovirus	Qualitative detection of Astrovirus in faeces in realtime	96 rxn
------------	------------------	-------------------	---	--------

Campylobacter

<i>ifp</i>	<i>PCAM0050</i>	Campylobacter	Qualitative detection of Campylobacter coli, laridis, jejuni in faeces in realtime (detection of the 16S rRNA gene; probe with FAM)	96 rxn
------------	-----------------	----------------------	---	--------

Clostridium difficile

<i>ifp</i>	<i>PCDI0090</i>	Clostridium difficile	Qualitative detection of Clostridium difficile toxin A and toxin B gene in faeces in realtime (probe with FAM reporter).	96 rxn
------------	-----------------	------------------------------	--	--------

Cryptosporidium / Giardia

<i>ifp</i>	<i>PCRY0120</i>	Cryptosporidium parvum	Qualitative detection of Cryptosporidium parvum in faeces in realtime (detection of the COPW gene; probe with FAM reporter).	96 rxn
------------	-----------------	-------------------------------	--	--------

Entamoeba

<i>ifp</i>	<i>PENT0100</i>	Entamoeba histolytica	Qualitative detection of Entamoeba histolytica in faeces (detection of the small subunit rRNA gene, probe with FAM reporter).	96 rxn
------------	-----------------	------------------------------	---	--------

Entero Hemorrhagic Escherichia coli

<i>ifp</i>	<i>PSTE0060</i>	VTEC stx1 / stx2	Qualitative detection of Shigatoxin producing Escherichia coli in faeces in realtime (detection of the stx1 and stx2 genes; probe with FAM reporter).	96 rxn
------------	-----------------	-------------------------	---	--------

Molecular Biology

PCR Fast Tests

Giardia				
<i>ifp</i>	<i>PGIA0080</i>	Giardia lamblia	Qualitative detection of Giardia lamblia in faeces in realtime (detection of the β -giardin gene; probe with FAM reporter).	96 rxn
Helicobacter pylori				
<i>ifp</i>	<i>PHEL0020</i>	Helicobacter pylori	Qualitative detection of Helicobacter pylori in faeces in realtime (detection of the urease C gene; probe with FAM reporter).	96 rxn
Listeria				
<i>ifp</i>	<i>PLIS0040</i>	Listeria monocytogenes	Qualitative detection of Listeria monocytogenes in faeces in realtime (detection of the prfA gene; probe with FAM reporter).	96 rxn
Microsporidia				
<i>ifp</i>	<i>PMSPO160</i>	Microsporidia	Qualitative detection of Microsporidia in faeces in realtime (detection of the SSU rRNA gene; probe with FAM reporter).	96 rxn
Norovirus				
<i>ifp</i>	<i>PNOR0110</i>	Norovirus	Qualitative detection of Norovirus genotype I and genotype II in faeces in realtime (probe with FAM reporter).	96 rxn
Rotavirus				
<i>ifp</i>	<i>PROT0200S</i>	Rotavirus	Qualitative detection of Rotavirus in faeces in realtime	96 rxn
Salmonella				
<i>ifp</i>	<i>PSAL0030</i>	Salmonella	Qualitative detection of Salmonella in faeces in realtime (detection of the invA gene; probe with FAM reporter).	96 rxn
Shigella				
<i>ifp</i>	<i>PSHI0070</i>	Shigella	Qualitative detection of Shigella in faeces in realtime (detection of the ipaH gene; probe with FAM reporter).	96 rxn
Yersinia enterocolitica				
<i>ifp</i>	<i>PYER0010</i>	Yersinia enterocolitica	Qualitative detection of Yersinia enterocolitica in faeces in realtime (detection of the ail gene; probe with FAM reporter).	96 rxn

Rapid Tests

Virology

Adenovirus				
Meridian Bioscience	751030 300651	ImmunoCard STAT! Adenovirus	One step immunochromatographic assay	30 tests
EBV				
Meridian Bioscience	751830 300041	TRU EBV-G	One step Immuno-Chromatographic assay for determination of anti-VCA IgG and anti-EBNA-1 IgG	30 tests
Meridian Bioscience	776500	Monolert	Rapid EIA test for anti-EBNA-1 IgG and IgM status	100 tests
Meridian Bioscience	776520	Monolert Control Sera	Positive and negative control sera	25 tests
Influenza Virus				
Meridian Bioscience	751230	TRU FLU	One step Immuno-Chromatographic assay, with separate A / B results	32 tests
Meridian Bioscience	751110	FLU/RSV Positive Control	3 vials of 8 tests each	24 tests
Norovirus				
Meridian Bioscience	757020 300752	ImmunoCard STAT! Norovirus	One step Immuno-Chromatographic Assay for the detection of Norovirus antigens in stool specimens	20 tests
Meridian Bioscience	757001 300759	ImmunoCard STAT! Norovirus positive control		1 test
Rotavirus				
Meridian Bioscience	750030 300551	ImmunoCard STAT! Rotavirus	One step immunochromatographic assay	30 tests
Meridian Bioscience	751120 300552	Rapid Strip Rota / Adeno	One step immunochromatographic assay with separate bands for Rotavirus and Adenovirus	20 tests
RSV				
Meridian Bioscience	751110	FLU/RSV Positive Control	3 vials of 8 tests each	24 tests
Meridian Bioscience	751330 300453	TRU RSV	One step Immuno-Chromatographic assay for the detection of RSV antigen	32 tests

Bacteriology

Campylobacter				
Meridian Bioscience	751530	ImmunoCard STAT! Campy	One step Immuno-Chromatographic Assay for the detection of Campylobacter antigens (C. jejuni, C. coli) in stool specimens	30 tests
Chlamydia				
Meridian Bioscience	757120	ImmunoCard STAT! Chlamydia pneumoniae IgM	Rapid one step test for Chlamydia pneumoniae IgM antibodies detection in serum or fingertip blood specimens	20 tests

Rapid Tests

Bacteriology

<i>Clostridium difficile</i>				
Meridian Bioscience	716050 303150	ImmunoCard C. difficile GDH	Rapid EIA for the detection of Clostridium difficile glutamate dehydrogenase (GDH) in fecal specimens	50 tests
Meridian Bioscience	712050 303151	ImmunoCard Toxins A&B	Rapid membrane EIA for the detection of C. diff. Toxins A&B	50 tests
<i>Entero Hemorrhagic Escherichia coli</i>				
Meridian Bioscience	750530 303202	ImmunoCard STAT! E. coli O157 PLUS	One step immunochromatographic assay test for the detection of E. coli O157 from stool specimens and from culture	30 tests
Meridian Bioscience	751630 303203	ImmunoCard STAT! EHEC	Lateral flow test for the detection of E. coli Shiga Toxins 1&2 (with differentiation) from enriched stool specimens	30 tests
Meridian Bioscience	7516MC2 303204	McConkey Broth	For use with ImmunoCard STAT! EHEC for enrichment of stool specimens	20 x 5ml
<i>Helicobacter pylori</i>				
Meridian Bioscience	750705 303050	ImmunoCard STAT! HpSA	One step immunochromatographic assay for the detection of H. pylori Specific Antigens in stool specimens	5 tests
Meridian Bioscience	750720 303052	ImmunoCard STAT! HpSA	One step immunochromatographic assay for the detection of H. pylori Specific Antigens in stool specimens	20 tests
Meridian Bioscience	710030 303055	ImmunoCard H. pylori	m-EIA for IgG detection in serum and whole blood	30 tests
<i>Legionella</i>				
Meridian Bioscience	751930	TRU Legionella	Rapid test for the qualitative detection of Legionella pneumophila serogroup 1 antigen in human urine specimens.	32 tests
<i>Mycoplasma pneumoniae</i>				
Meridian Bioscience	709030 303551	ImmunoCard Mycoplasma	Rapid membrane EIA for M. pneumoniae IgM detection	30 tests
<i>Streptococcus</i>				
Teomed	241010	Teomed Strep A Test	Rapid detection of Group A Streptococcal antigen from patient throat swab specimens using flocced swabs	10 tests
Teomed	241020	Teomed Strep A Test	Rapid detection of Group A Streptococcal antigen from patient throat swab specimens using flocced swabs	20 tests

Parasitology

<i>Cryptosporidium / Giardia</i>				
Meridian Bioscience	750830 306151	ImmunoCard STAT! Crypto / Giardia	Rapid immunoassay for the qualitative detection of Cryptosporidium parvum and Giardia lamblia specific antigens from stool specimens	30 tests
Meridian Bioscience	750820	ImmunoCard STAT! Crypto / Giardia Positive Controls	1 ml Giardia and 1 ml Cryptosporidium Positive Control vials	16 tests

Rapid Tests

Parasitology

Cryptosporidium / Giardia

<i>Meridian Bioscience</i>	757220	Rapid Strip Crypto & Giardia	One step Immuno-Chromatographic Assay for the detection of <i>Cryptosporidium parvum</i> and <i>Giardia lamblia</i> antigens in stool	20 tests
--------------------------------	--------	---	---	----------

Clinical Chemistry

Occult Blood in Stool

<i>Epitope Diagnostics</i>	KT313 230010	EpiTuub Fecal Occult Blood Test	Monoclonal based immunochromatographic assay for the detection of fecal occult blood	30 tests
--------------------------------	-----------------	--	--	----------

ELISA - Antibody Detection

Virology

Adenovirus				
NovaTec	ADVA0010	Adenovirus IgA	Highly purified antigens of strain Adenoid 6	96 tests
NovaTec	ADVG0010	Adenovirus IgG	Highly purified antigens of strain Adenoid 6	96 tests
NovaTec	ADVM0010	Adenovirus IgM	Highly purified antigens of strain Adenoid 6	96 tests
Chikungunya Virus				
NovaTec	CHIG0590	Chikungunya IgG	Chikungunya antigens	96 tests
NovaTec	CHIM0590	Chikungunya IgM	μ-capture, Chikungunya antigens	96 tests
CMV				
NovaTec	CMVG0110	Cytomegalovirus (CMV) IgG	CMV strain AD-169 antigens	96 tests
NovaTec	CMVM0110	Cytomegalovirus (CMV) IgM	CMV strain Davis antigens	96 tests
Zeus Scientific	9Z9501G	CMV IgG	Detection of IgG class antibodies to cytomegalovirus (Cytomegalovirus) in human serum	96 tests
Zeus Scientific	9Z9501M	CMV IgM	Detection of IgM class antibodies to cytomegalovirus (Cytomegalovirus) in human serum	96 tests
Dengue Virus				
NovaTec	DENG0120	Dengue IgG	Purified Dengue Virus 2 antigens from strain 16681	96 tests
NovaTec	DENM0120	Dengue IgM	Purified Dengue Virus 2 antigens from strain 16681	96 tests
EBV				
Zeus Scientific	9Z9201G 360001	EBV-VCA IgG EIA	Based on affinity purified EBV VCA gp125 antigen	96 tests
Zeus Scientific	9Z9201M 360002	EBV-VCA IgM EIA	Based on affinity purified EBV VCA gp125 antigen	96 tests
Zeus Scientific	9Z9461G 360003	EBV-NA IgG EIA	Based on synthetic peptide containing active sites specific for the detection of IgG antibodies to EBNA-1	96 tests
Zeus Scientific	9Z9481G 360004	EBV-EA IgG EIA	Based on semi-purified EBV-EA antigen derived from induced Raji cells	96 tests
FSME (TBE; Tick-Borne Encephalitis)				
Progen Biotechnik	7701010 202001	IMMUNOZYM FSME (TBE) IgG	For the quantitative detection of IgG antibodies against TBE virus in serum, plasma and cerebrospinal fluid	96 tests
Progen Biotechnik	7701045 202002	IMMUNOZYM FSME (TBE) IgM	For the quantitative detection of IgM antibodies against TBE virus in serum, plasma and cerebrospinal fluid	96 tests
Progen Biotechnik	7701075	IMMUNOZYM FSME (TBE) IgG All Species	For the quantitative detection of antibodies against TBE virus in different species, e.g. dog, horse, cattle.	96 tests

ELISA - Antibody Detection

Virology

Hanta Virus				
Progen Biotechnik	PR59165	Hantavirus (Dobrava/Hantaan) IgG/IgM	Semiquantitative Elisa for the detection of IgG and IgM antibodies against hantavirus (Dobrava and Hantaan serotypes) in human serum	96 tests
Progen Biotechnik	PR59156	Hantavirus (Puumala) IgG/IgM	Semiquantitative Elisa for the detection of IgG and IgM antibodies against hantavirus (Puumala group) in human serum	96 tests
HSV				
NovaTec	HSV1G0500 320704	Herpes simplex virus Type 1 (HSV-1) IgG	recombinant gG1 protein of HSV Type 1	96 tests
NovaTec	HSV1M0500 320705	Herpes simplex virus Type 1 (HSV-1) IgM	recombinant gG1 protein of HSV Type 1	96 tests
NovaTec	HSV2G0540 320714	Herpes simplex virus Type 2 (HSV-2) IgG	recombinant gG2 protein of HSV Type 2	96 tests
NovaTec	HSV2M0540 320715	Herpes simplex virus Type 2 (HSV-2) IgM	recombinant gG2 protein of HSV Type 2	96 tests
NovaTec	HSVG0250	Herpes simplex virus Type 1+2 (HSV-1+2) IgG	purified membrane antigens of strain Mac Intyre (ATCC VR-539) (HSV1) and strain Ms (ATCC VR-540) (HSV2)	96 tests
NovaTec	HSV M0250 320725	Herpes simplex virus Type 1+2 (HSV-1+2) IgM	purified membrane antigens of strain Mac Intyre (ATCC VR-539) (HSV1) and strain Ms (ATCC VR-540) (HSV2)	96 tests
Zeus Scientific	9Z9761G	HSV gG-1 IgG	Type specific detection of IgG antibodies to HSV-1	96 tests
Zeus Scientific	9Z9781G	HSV gG-2 IgG	Type specific detection of IgG antibodies to HSV-2	96 tests
Zeus Scientific	9Z9771M	HSV-1&2 IgM	Detection of IgM class antibodies to Herpes Simplex Virus (HSV)	96 tests
Influenza Virus				
NovaTec	INFA0290	Influenza virus A IgA	Purified Influenza A grade 2 antigens (strain Texas 1/77 (H3N2))	96 tests
NovaTec	INFG0290	Influenza virus A IgG	Purified Influenza A grade 2 antigens (strain Texas 1/77 (H3N2))	96 tests
NovaTec	INFM0290	Influenza virus A IgM	Purified Influenza A grade 2 antigens (strain Texas 1/77 (H3N2))	96 tests
NovaTec	INFA0300	Influenza virus B IgA	Purified Influenza B grade 2 antigens (strain Hong Kong 5/72)	96 tests
NovaTec	INFG0300	Influenza virus B IgG	Purified Influenza B grade 2 antigens (strain Hong Kong 5/72)	96 tests
NovaTec	INFM0300	Influenza virus B IgM	Purified Influenza B grade 2 antigens (strain Hong Kong 5/72)	96 tests

ELISA - Antibody Detection

Virology

Measles (Rubeola) Virus				
NovaTec	MEAG0330	Measles virus IgG	Purified measles virus antigens of strain Edmonton (ATCC VR-24)	96 tests
NovaTec	MEAM0330 321212	Measles virus IgM	Purified measles virus antigens of strain Edmonton (ATCC VR-24)	96 tests
Zeus Scientific	9Z9271G	Measles (Rubeola) IgG	Detection of IgG antibodies to Measles (Rubeola) virus	96 tests
Mumps Virus				
NovaTec	MUMG0340 321201	Mumps virus IgG	Purified Parotitis Virus antigens, strain Enders (ATCC VR-106)	96 tests
NovaTec	MUMM0340 321202	Mumps virus IgM	Purified Parotitis Virus antigens, strain Enders (ATCC VR-106)	96 tests
Zeus Scientific	9Z9281G	Mumps IgG	Detection of Mumps IgG antibody	96 tests
Parainfluenza Virus				
NovaTec	PAIA0360	Parainfluenza 1,2,3 IgA	Purified Parainfluenza virus antigens	96 tests
NovaTec	PAIG0360	Parainfluenza 1,2,3 IgG	Purified Parainfluenza virus antigens	96 tests
Parvovirus				
NovaTec	PARG0370 321511	Parvovirus B19 IgG	recombinant Parvovirus B19 antigens	96 tests
NovaTec	PARM0370 321512	Parvovirus B19 IgM	recombinant Parvovirus B19 antigens	96 tests
RSV				
NovaTec	RSVA0380	Respiratory syncytial virus (RSV) IgA	Purified Respiratory syncytial Virus antigens (strain RuS)	96 tests
NovaTec	RSVG0380	Respiratory syncytial virus (RSV) IgG	Purified Respiratory syncytial Virus antigens (strain RuS)	96 tests
NovaTec	RSVM0380	Respiratory syncytial virus (RSV) IgM	Purified Respiratory syncytial Virus antigens (strain RuS)	96 tests
Rubella Virus				
NovaTec	RUBG0400	Rubella virus IgG	quantitative, Rubella K1S antigens	96 tests
NovaTec	RUBM0400	Rubella virus IgM	μ-capture, purified rubella antigens	96 tests
Zeus Scientific	9Z9801G	Rubella IgG	Detection of IgG antibodies to rubella virus	96 tests
Zeus Scientific	9Z9801M	Rubella IgM	Detection of IgM antibodies to rubella virus	96 tests
VZV				
NovaTec	VZVA0490 322103	Varizella-Zoster (VZV) IgA	Purified Varicella-Zoster-Virus antigens	96 tests

ELISA - Antibody Detection

Virology

VZV				
NovaTec	VZVG0490 322101	Varizella-Zoster (VZV) IgG	Purified Varicella-Zoster-Virus antigens	96 tests
NovaTec	VZVM0490 322102	Varizella-Zoster (VZV) IgM	Purified Varicella-Zoster-Virus antigens	96 tests
Zeus Scientific	9Z9331G	Varicella-Zoster (VZV) IgG	Detection of Varicella Zoster Virus (VZV) IgG antibodies	96 tests
Zeus Scientific	9Z9331M	Varicella-Zoster (VZV) IgM	Detection of Varicella Zoster Virus (VZV) IgM antibodies	96 tests

Bacteriology

<i>Bordetella pertussis</i>				
NovaTec	BOPA0030 324101	Bordetella pertussis IgA	B. pertussis Toxin and filamentous Hemagglutinine antigens	96 tests
NovaTec	BOPG0030 324102	Bordetella pertussis IgG	B. pertussis Toxin and filamentous Hemagglutinine antigens	96 tests
NovaTec	BOPM0030 324103	Bordetella pertussis IgM	Purified lysate antigens of B. pertussis	96 tests

<i>Borrelia</i>				
NovaTec	BORG0040	Borrelia burgdorferi IgG (Lyme)	Recombinant antigens of B. garinii and B. afzelii (OspC, p100, p18, p41i)	96 tests
NovaTec	BORM0040	Borrelia burgdorferi IgM (Lyme)	Recombinant antigens of B. garinii and B. afzelii (OspC, p41i)	96 tests

<i>Brucella</i>				
NovaTec	BRUG0050	Brucella IgG	Purified Brucella abortus (W 99) antigens	96 tests
NovaTec	BRUM0050	Brucella IgM	Purified Brucella abortus (W 99) antigens	96 tests

<i>Chlamydia</i>				
NovaTec	CHLA0070	Chlamydia trachomatis IgA	Highly purified C. trachomatis LGV Type II strain 434 antigens	96 tests
NovaTec	CHLG0070	Chlamydia trachomatis IgG	Highly purified C. trachomatis LGV Type II strain 434 antigens	96 tests
NovaTec	CHLM0070	Chlamydia trachomatis IgM	Highly purified C. trachomatis LGV Type II strain 434 antigens	96 tests
NovaTec	CHLA0510	Chlamydia pneumoniae IgA	Highly purified C. pneumoniae CWL-029 antigens	96 tests
NovaTec	CHLG0510	Chlamydia pneumoniae IgG	Highly purified C. pneumoniae CWL-029 antigens	96 tests
NovaTec	CHLM0510	Chlamydia pneumoniae IgM	Highly purified C. pneumoniae CWL-029 antigens	96 tests

ELISA - Antibody Detection

Bacteriology

<i>Clostridium tetani</i>				
Progen Biotechnik	7702010 202003	IMMUNOZYM Tetanus	Quantitative determination of anti-tetanus toxoid antibodies in human serum and human plasma. Suitable for control of humoral immune status after vaccination ("Vaccination Management")	96 tests
<i>Corynebacterium diphtheriae</i>				
Progen Biotechnik	7702040 202004	IMMUNOZYM Diphtheria	Quantitative determination of anti-diphtheria toxoid antibodies in human serum and human plasma. Suitable for control of humoral immune status after vaccination ("Vaccination Management")	96 tests
<i>Coxiella burnetii</i>				
NovaTec	COX2G0600	Coxiella Burnetii (Q-Fever) Phase 1 IgG	Qualitative determination of IgG antibodies against Coxiella burnetii phase 1 in human serum	96 tests
NovaTec	COX1G0600	Coxiella Burnetii (Q-Fever) Phase 2 IgG	Qualitative determination of IgG antibodies against Coxiella burnetii phase 2 in human serum	96 tests
NovaTec	COX2M0600	Coxiella Burnetii (Q-Fever) Phase 2 IgM	Qualitative determination of IgM antibodies against Coxiella burnetii phase 2 in human serum	96 tests
NovaTec	COX1G0600	Coxiella Burnetii (Q-Fever) Phase 1 IgG	Qualitative determination of IgG antibodies against Coxiella burnetii phase 1 in the early stages of infection in human serum	96 tests
<i>Haemophilus influenzae type B (HiB)</i>				
Progen Biotechnik	7704010 202005	IMMUNOZYM HiB IgG	Quantitative determination of IgG antibodies against Haemophilus influenzae type B in human serum and human plasma. Suitable for control of humoral immune status after vaccination ("Vaccination Management")	96 tests
<i>Helicobacter pylori</i>				
NovaTec	HELA0220	Helicobacter pylori IgA	Quantitative, recombinant CagA- and recombinant urease-antigens	96 tests
NovaTec	HELG0220 324701	Helicobacter pylori IgG	Quantitative, highly purified proteins associated with CagA genes (120 KD) and VaCA genes (87 KD) as well as urease-antigens	96 tests
Zeus Scientific	3Z51021G	Helicobacter pylori IgG	Detection of IgG antibodies to Helicobacter pylori	96 tests
<i>Legionella</i>				
Meridian Bioscience	633096	Premier L. pneumophila	EIA for the detection of antibodies against L. pneumophila (serogr. 1-6) in human serum	96 tests
Zeus Scientific	3Z15051	Legionella IgG/M/A	Qualitative detection of total antibody (IgG/IgM/IgA) to Legionella pneumophila serogroups 1-6 in human sera.	96 tests
<i>Mycobacteria</i>				
Progen Biotechnik	311905	AndaELISA Mycobacteria IgG	For the quantitative determination of IgG antibodies against mycobacterium tuberculosis and other mycobacteria in human serum	96 tests

ELISA - Antibody Detection

Bacteriology

Mycobacteria				
Progen Biotechnik	311907	AndaELISA Mycobacteria IgA	For the quantitative determination of IgA antibodies against mycobacterium tuberculosis and other mycobacteria in human serum	96 tests
Progen Biotechnik	311915	AndaELISA Mycobacteria IgM	For the quantitative determination of IgM antibodies against mycobacterium tuberculosis and other mycobacteria in human serum	96 tests
Mycoplasma pneumoniae				
NovaTec	MYCA0350	Mycoplasma pneumoniae IgA	Purified Mycoplasma pneumoniae antigens (strain FN)	96 tests
Zeus Scientific	3Z17601G	Mycoplasma IgG	Detection of IgG antibodies to Mycoplasma pneumoniae (strain FN)	96 tests
Zeus Scientific	3Z17601M	Mycoplasma IgM	Detection of IgM antibodies to Mycoplasma pneumoniae (strain FN)	96 tests
Syphilis				
NovaTec	TREG0470	Treponema pallidum Ab	Recombinant antigens	96 tests
Zeus Scientific	3Z7601	T. pallidum (Syphilis) IgG	Detection of IgG class antibodies to 17 kDa antigen from T. pallidum	96 tests

Parasitology

Ascaris				
NovaTec	ASCG0020	Ascaris lumbricoides IgG	PBS extract of Ascaris	96 tests
Leishmania				
NovaTec	LEIG0310	Leishmania IgG	Purified Leishmania antigens	96 tests
NovaTec	LEIG0310V	Leishmania IgG canine	Purified Leishmania antigens	96 tests
Toxoplasma gondii				
NovaTec	TOXM0460	Toxoplasma gondii IgM	μ-Capture, P30 antigen	96 tests
NovaTec	TOXG0460	Toxoplasma gondii IgG	quantitative, Toxoplasma gondii (Lysate) antigens	96 tests
NovaTec	TOXGA460	Toxoplasma gondii IgG Avidity Reagents	IgG avidity reagents intended to differentiate between past and acute infection. To be used with TOXG0460 Toxoplasma gondii IgG	96 tests
Zeus Scientific	8Z8651G	Toxoplasma IgG	Detection of IgG antibodies to Toxoplasma gondii	96 tests
Zeus Scientific	8Z8651M	Toxoplasma IgM	Detection of IgM antibodies to Toxoplasma gondii	96 tests

ELISA - Antibody Detection

Parasitology

<i>Trichinella</i>				
NovaTec	TRIG0480	Trichinella spiralis IgG	Trichinella antigens	96 tests
<i>Trypanosoma Cruzi (Chagas)</i>				
NovaTec	CHAG0560	Chagas (Trypanosoma cruzi) IgG	Recombinant Trypanosoma cruzi antigens (TcF)	96 tests

Mycology

<i>Candida albicans</i>				
NovaTec	CANA0060	Candida albicans IgA	Purified antigens from Candida albicans	96 tests
NovaTec	CANG0060	Candida albicans IgG	Purified antigens from Candida albicans	96 tests
NovaTec	CANM0060	Candida albicans IgM	Purified antigens from Candida albicans	96 tests
<i>Fungal Serology</i>				
Meridian Bioscience	603096	Premier Coccidioides	Detection of Anti-C. immitis IgG and IgM	96 tests

Immunology

<i>Autoimmune Diagnostics</i>				
Zeus Scientific	2Z29001G	ANA Screen		96 tests
Zeus Scientific	2Z2851G	ENA Screen		96 tests
Zeus Scientific	2Z2801G	ENA Profile-6, semi-quant.		96 tests
Zeus Scientific	2Z2881G	dsDNA, semi-quant.		96 tests
Zeus Scientific	2Z2811G	SSA (Ro), semi-quant.		96 tests
Zeus Scientific	2Z2821G	SSB (La), semi-quant.		96 tests
Zeus Scientific	2Z2831G	Sm, semi-quant.		96 tests
Zeus Scientific	2Z2841G	Sm/RNP, semi-quant.		96 tests
Zeus Scientific	2Z2871G	Jo-1, semi-quant.		96 tests
Zeus Scientific	2Z2861G	Scl-70, semi-quant.		96 tests
Zeus Scientific	2Z9681G	ANCA Screen (PR3/MPO)		96 tests

ELISA - Antibody Detection

Immunology

<i>Autoimmune Diagnostics</i>			
Zeus Scientific	2Z9671G	pANCA (MPO), semi-quant.	96 tests
Zeus Scientific	2Z9691G	cANCA (PR3), semi-quant.	96 tests
Zeus Scientific	2Z51051	Cardiolipin Screen (IgG/IgM), semi-quant.	96 tests
Zeus Scientific	2Z51051A	Cardiolipin IgA, semi-quant.	96 tests
Zeus Scientific	2Z51051M	Cardiolipin IgM, semi-quant.	96 tests
Zeus Scientific	2Z51051G	Cardiolipin IgG, semi-quant.	96 tests
Zeus Scientific	2Z61001A	Gliadin IgA, semi-quant.	96 tests
Zeus Scientific	2Z61001G	Gliadin IgG, semi-quant.	96 tests
Zeus Scientific	2Z5061G	Thyroglobulin (Tg), semi-quant.	96 tests
Zeus Scientific	2Z5051G	Thyroperoxidase (TPO), semi-quant.	96 tests
Zeus Scientific	2Z9291M	Rheumatoid Factor (RF) IgM, semi-quant.	96 tests

ELISA - Antigen Detection

Virology

Adenovirus				
Meridian Bioscience	696006	Premier Adenoclone 40/41	Microtiter EIA for direct detection of Adenoviruses serotypes 40/41 in stool samples	48 tests
Meridian Bioscience	696007	Premier Adenoclone	Microtiter EIA for detection of adenovirus antigen in stool samples and for confirmation in cell culture	48 tests
Rotavirus				
Meridian Bioscience	696004	Premier Rotaclone	Using Rotavirus monoclonal antibodies	48 tests

Bacteriology

Campylobacter				
Meridian Bioscience	618096	Premier Campy	Detection of Campylobacter antigens (C. jejuni, C. Coli) in stool specimens	96 tests
Clostridium difficile				
Meridian Bioscience	611096 303100	Premier C. difficile GDH	Detection of C. difficile Common Antigen in stool specimens	96 tests
Meridian Bioscience	616096 303101	Premier Toxins A&B	Microtiter EIA for the detection of C. diff. Toxins A&B in stool specimens	96 tests
Entero Hemorrhagic Escherichia coli				
Meridian Bioscience	608096 303201	Premier EHEC	Microtiter EIA for the detection of Shiga-like Toxins I & II in stool and food specimens	96 tests
Helicobacter pylori				
Meridian Bioscience	601396 303002	Premier Platinum HpSA PLUS	Monoclonal based microtiter EIA for the detection of H. pylori Specific Antigens in stool specimens	96 tests

Parasitology

Cryptosporidium / Giardia				
NovaTec	GIA0160S	Giardia lamblia Antigen in stool		96 tests

Mycology

Cryptococcus neoformans				
Meridian Bioscience	602096	Premier Cryptococcal Antigen	Microtiter EIA for the detection of Cr.neoformans antigen in serum and CSF	192 tests

ELISA - Antigen Detection

Clinical Chemistry

Lipid Metabolism

Progen Biotechnik	7700040 202004	IMMUNOZYM Lp(a)	Quantitative Elisa to determine the concentration of Lp(a)	96 tests
-------------------	-------------------	------------------------	--	----------

Occult Blood in Stool

Epitope Diagnostics	KT-850	Fecal Human Hemoglobin		96 tests
------------------------	--------	-------------------------------	--	----------

Immunoblot

Virology

Rubella Virus

NovaTec	RUBG2400	Rubella Virus IgG	Immunoblot for the qualitative detection of rubella virus specific IgG antibodies in human serum as a tool to discriminate between acute and past infections	10 tests
---------	----------	--------------------------	--	----------

Immunology

Allergy Diagnostics

Mediwiss Analytic	120210	AllergyScreen INH 16+T	Immunoblot for the determination of 16 different specific IgE to inhalative allergens (Swiss Panel)	12 tests
Mediwiss Analytic	120220	AllergyScreen NM 16+T	Immunoblot for the determination of 16 different specific IgE to alimentary allergens (Swiss Panel)	12 tests
Mediwiss Analytic	120230	AllergyScreen Panel INH 10	Immunoblot for the determination of 10 different specific IgE to inhalative allergens	12 tests
Mediwiss Analytic	120400	AllergyScreen Panel NM/INH 10	Immunoblot for the determination of specific IgE to 7 alimentary and 3 inhalative allergens	12 tests
Mediwiss Analytic	122005	AllergyScreen Panel 1	Immunoblot for the determination of 20 different specific IgE to inhalative and alimentary allergens	12 tests
Mediwiss Analytic	122010	AllergyScreen Panel 2	Immunoblot for the determination of 20 different specific IgE to inhalative allergens	12 tests
Mediwiss Analytic	122020	AllergyScreen Panel 3	Immunoblot for the determination of 20 different specific IgE to alimentary allergens	12 tests
Mediwiss Analytic	122030	AllergyScreen Panel 4	Immunoblot for the determination of 20 different specific IgE to inhalative and alimentary allergens. Designed for pediatric patients.	12 tests
Mediwiss Analytic	122040	AllergyScreen Panel 5	Immunoblot for the determination of 20 different specific IgE to inhalative allergens and bee/wasp venom.	12 tests
Mediwiss Analytic	120401	Improvio C Reader	Scanner and software for the analysis of the AllergyScreen tests	1 pce tests

IFA - Antibody Detection

Virology

CMV				
Zeus Scientific	FA9001G	CMV IgG	Qualitative detection of IgG class antibodies to cytomegalovirus (CMV) in human serum	10x10 tests
EBV				
Meridian Bioscience	EB-100 300011	Merifluor EBV IgG	Slides coated with P3-HR1 cells	10x10 tests
Meridian Bioscience	EB-150 300012	Merifluor EBV IgM	Slides coated with P3-HR1 cells	10x10 tests
Meridian Bioscience	EN-100 300013	Merifluor EBNA ACIF	Slides coated with 20-40% EBNA pos. Burkitt lymphoma Raji cell line and with EBNA neg. cells; conjugate contains Evans Blue	10x10 tests
Meridian Bioscience	EA-101 300014	Merifluor EBV-EA	Slides coated with Raji cells from Burkitt lymphoma	10x10 tests
Meridian Bioscience	EB-211 300021	EBV IgM/IgG slides pack	Slides coated with P3-HR1 cells	10x10 tests
Meridian Bioscience	EA-211 300024	EBV-EA slides pack	Slides coated with Raji cells from Burkitt lymphoma	10x10 tests
Hanta Virus				
Progen Biotechnik	PR77065 202504	Hantavirus (Dobrava/Hantaan) Antibody IFA	For the detection of antibodies against Hantavirus Hantaan and Dobrava serotypes in human serum.	60 tests
Progen Biotechnik	PR77056 202505	Hantavirus (Puumala) Antibody IFA	For the detection of antibodies against hantavirus (Puumala group) in human serum	60 tests
Progen Biotechnik	PR77060 202506	Hantavirus (Seoul) Antibody IFA	For the detection of antibodies against hantavirus (Seoul group) in human serum	60 tests
Progen Biotechnik	FR77090 203110	Hantavirus IF Reagent Kit Rat/Mouse	For use in an indirect immunofluorescence test for the the determination of rat or mouse antibodies against hantavirus (hantaan/puumala/seoul) in rat and mouse serum.	60 tests
Progen Biotechnik	PR97065	Slide Pack (Hantaan)		60 tests
Progen Biotechnik	PR97056	Slide pack (Puumala)		60 tests
Progen Biotechnik	PR97060 203111	Slide pack (Seoul)		60 tests
HSV				
Meridian Bioscience	HS-100 300111	Merifluor HSV IgG	Slides coated with uninfected control cells and HSV infected cells	10x10 tests
Meridian Bioscience	HS-150 300112	Merifluor HSV IgM	Slides coated with uninfected control cells and HSV infected cells	10x10 tests

IFA - Antibody Detection

Virology

Measles (Rubeola) Virus				
Zeus Scientific	FA9251G	Measles (Rubeola) IgG	Quantitative and qualitative detection of Measles IgG antibodies	10x10 tests
VZV				
Meridian Bioscience	VZ-150 300212	Merifluor VZV IgM	Slides coated with VZV (Ellen strain) infected human fibroblast cells and uninfected control cells	10x10 tests
Meridian Bioscience	VZ-100 300211	Merifluor VZV IgG	Slides coated with VZV (Ellen strain) infected human fibroblast cells and uninfected control cells	10x10 tests

Bacteriology

Legionella				
Meridian Bioscience	4873171 303411	L. pneumophila serogroups 1-8	IFA for the detection of antibodies to L.pneumophila serogroups 1-8, 20 slides (wells 1-5 are coated with inactivated bacteria of serogroup 1, 4, 6, 8; wells 6-10 with inact. bacteria of serogroup 2, 3, 5, 7)	100 tests
Meridian Bioscience	4884161	L. pneumophila polyvalent Serogroups 1-6		10 tests
Meridian Bioscience	4873891	L. pneumophila polyvalent Serogroups 1-8	(Pool I (1, 4, 6, 8), Pool II (2, 3, 5, 7))	10 tests
Meridian Bioscience	4873001	L. pneumophila serogroup 1		10 tests
Meridian Bioscience	4873002	L. pneumophila serogroup 2		10 tests
Meridian Bioscience	4873003	L. pneumophila serogroup 3		10 tests
Meridian Bioscience	4873004	L. pneumophila serogroup 4		10 tests
Meridian Bioscience	4873005	L. pneumophila serogroup 5		10 tests
Meridian Bioscience	4873006	L. pneumophila serogroup 6		10 tests
Meridian Bioscience	4873141	L.pneumophila negative control serum		1 ml
Meridian Bioscience	4873011	L.pneumophila positive control serum		1 ml
Mycoplasma pneumoniae				
Zeus Scientific	FA17002G-CT	Mycoplasma pneumonia IgG Crowntitre	Qualitative and semi-quantitative detection of Mycoplasma pneumoniae IgG antibodies	10x10 tests
Zeus Scientific	FA17002M-CT	Mycoplasma pneumonia IgM Crowntitre	Qualitative and semi-quantitative detection of Mycoplasma pneumoniae IgM antibodies	10x10 tests

IFA - Antibody Detection

Bacteriology

Syphilis				
Zeus Scientific	FA7001	FTA-ABS (Syphilis) T. pallidum	Qualitative determination of antibodies to Treponema pallidum	10x10 tests
Zeus Scientific	FA7001DS	FTA-ABS (Syphilis) T. pallidum Double Stain	Designed to confirm positive non-treponemal reagin tests for syphilis	10x10 tests
Zeus Scientific	FA7001T	FTA-ABS (Syphilis) T. pallidum Titratable	Intended to be used as an aid in the confirmation of syphilis antibodies	10x10 tests

Parasitology

Toxoplasma gondii				
Zeus Scientific	FA8001G	Toxoplasma IgG	Detection of Toxoplasma gondii IgG antibodies	120 Tests
Zeus Scientific	FA8001M	Toxoplasma IgM	Detection of Toxoplasma gondii IgM antibodies	120 Tests

Immunology

Autoimmune Diagnostics				
Zeus Scientific	2401	HEp-2 IFA liquid	Prestandardized kit utilizing tissue cell culture for the qualitative and semi-quantitative detection of anti-nuclear antibodies (ANA)	10x5 wells
Zeus Scientific	1001	nDNA C. luciliae	Qualitative and semi-quantitative detection of antibodies to native DANN	6x8 wells
Zeus Scientific	2001	ANA Rat Liver Tissue	Qualitative and semi-quantitative detection of antinuclear antibodies	9x6 wells
Zeus Scientific	3001	AMA Rat Kidney Tissue	Qualitative and semi-quantitative detection of mitochondrial antibodies	8x6 wells
Zeus Scientific	2201-8	ANA Mouse Kidney Tissue	Qualitative and semi-quantitative detection of anti-nuclear antibodies	7x8 wells
Zeus Scientific	6001	ASA Monkey Esophagus Tissue	Qualitative and semi-quantitative detection of antibodies associated with pemphigus and bullous pemphigoid	8x6 wells
Zeus Scientific	4001	ASMA Rat Stomach Tissue	Qualitative and semi-quantitative detection of smooth muscle antibodies	8x6 wells
Zeus Scientific	5001	TA Monkey Thyroid Tissue	Qualitative and semi-quantitative detection of thyroglobulin antibodies	8x6 wells
Zeus Scientific	3421	Autoantibody Screen Mouse Kidney / Stomach / Liver	Qualitative and semi-quantitative detection of antinuclear, mitochondrial, smooth muscle & parietal antibodies	6x8 wells

IFA - Antibody Detection

Immunology

Autoimmune Diagnostics

Zeus Scientific	3401-SP	Autoantibody Screen Rat Kidney / Stomach / Liver	Qualitative and semi-quantitative detection of antinuclear, smooth muscle & parietal cell antibodies	6x8 wells
Zeus Scientific	3401	Autoantibody Screen Rat Kidney / Stomach	Qualitative and semi-quantitative detection of antinuclear, mitochondrial, smooth muscle & parietal antibodies	6x8 wells

IFA - Antigen Detection

Virology

VZV

Meridian	500113	Merifluor VZV	For detection of VZV in specimens and culture	60 tests
Bioscience	300210			

Bacteriology

Legionella

Zeus Scientific	12001-6	Legionella DFA	Direct Fluorescent Antibody detection reagents are designed for the detection of Legionella organisms in various patient specimens	1 set
-----------------	---------	-----------------------	--	-------

Parasitology

Cryptosporidium / Giardia

Meridian	250050	Merifluor Crypto & Giardia	Direct FA for Cryptosporidium and Giardia in stool specimens	50 tests
Bioscience	306110			

Mycology

Pneumocystis carinii

Meridian	222030	Merifluor Pneumocystis	Antigen detection from respiratory specimens	30 tests
Bioscience	308211			

Agglutination Tests

Virology

		EBV		
Meridian	776060	Monospot	Slide test using fresh horse erythrocytes	20 tests
Bioscience	300051			
Meridian	776150	Monospot Latex	Latex particle agglutination assay	50 tests
Bioscience	300052			

Bacteriology

		<i>Yersinia enterocolitica</i>		
Progen Biotechnik	PR87025 202508	Yersinia enterocolitica 0:3	Monoclonal antibody for identification of bacteria colonies of <i>Yersinia enterocolitica</i> serogroup 0:3	40 tests
Progen Biotechnik	PR87035 202509	Yersinia enterocolitica 0:9	Monoclonal antibody for identification of bacteria colonies of <i>Yersinia enterocolitica</i> serogroup 0:9	40 tests

Mycology

		<i>Cryptococcus neoformans</i>		
Meridian	140100	Cryptococcal Latex Aggl. System (CALAS)	For the detection of <i>Cr. neoformans</i> Ag in serum and CSF	100 tests
Bioscience	308151			
Meridian	140050	CALAS Pronase Reagent	Enzyme for the treatment of serum samples	5 vials
Bioscience				

Collection and Processing Systems for Stool Specimens

Closed Collection & Concentration Systems

		Parapak		
Meridian Bioscience	970120	Para-Pak Macro-CON	stool concentration system	120 units
Meridian Bioscience	900312 307001	Para-Pak Clean vial	empty vial	120 units
Meridian Bioscience	900412	Para-Pak 10% Formalin	10% neutral buffered formalin	120 units
Meridian Bioscience	900212 307002	Para-Pak SAF	sodium acetate formalin	120 units
Meridian Bioscience	901312	Para-Pak Ecofix	new safer fixative without formalin	120 units
Meridian Bioscience	900112	Para-Pak LV-PVA	polyvinyl alcohol with mercuric chloride	120 units
Meridian Bioscience	901012	Para-Pak Zn-PVA	mercury free PVA (with zinc sulphate)	120 units
Meridian Bioscience	900612 307003	Para-Pak C&S	transport medium for enteric pathogens	120 units

Collection Systems with Integrated Filters

		Parapak Plus+Ultra		
Meridian Bioscience	950205	Para-Pak PLUS SAF	minimal specimen handling system containing fixative	50 units
Meridian Bioscience	950405	Para-Pak PLUS 10% Formalin	minimal specimen handling system containing fixative	50 units
Meridian Bioscience	950420	Para-Pak PLUS 10% Formalin	minimal specimen handling system containing fixative	200 units
Meridian Bioscience	950605	Para-Pak PLUS Ecofix	minimal specimen handling containing non-toxic fixative	50 units
Meridian Bioscience	950620	Para-Pak PLUS Ecofix	minimal specimen handling containing non-toxic fixative	200 units
Meridian Bioscience	980212 307021	Para-Pak ULTRA SAF	closed system for 15/50 ml tubes	120 units
Meridian Bioscience	980412	Para-Pak ULTRA 10% Formalin	closed system for 15/50 ml tubes	120 units
Meridian Bioscience	981312	Para-Pak ULTRA Ecofix	closed system for 15/50 ml tubes, safe fixative	120 units
Meridian Bioscience	800400	Bulk tubes and caps 15 ml	recommended for use with Para-Pak ULTRA	500 units

Filtering Centrifuge Tube

		SpinCon		
Meridian Bioscience	972200 307031	SpinCon	Filtering tube for easy concentration, no chemical required	200 units
Meridian Bioscience	972500 307032	SpinCon	Filtering tube for easy concentration, no chemical required	500 units

Collection and Processing Systems for Stool Specimens

Reagents for Parasitology

		Reagents		
Meridian Bioscience	401301	Ecofix bulk	new safer fixative without formalin	500 ml
Meridian Bioscience	FK4318 307044	Ethyl Acetate	recommended for the Ritchie concentration	225 ml
Meridian Bioscience	801300 307043	Hemo-De	safe substitute for Ethyl Acetate and Xylene	400 ml
Meridian Bioscience	400101 307041	Wheatley's Trichrome Stain	ready-to-use reagent for permanent staining	500 ml
Meridian Bioscience	400501 307042	Para-Pak Trichrome Blue Stain	specific, ready-to-use stain for Microsporidia	500 ml
Meridian Bioscience	801400	Ecostain Set	reagents set for staining Ecofix preparations	1 set

Other

Mycology

		Fungal Serology		
Meridian Bioscience	100201 308031	Histoplasma Antigen		1 ml
Meridian Bioscience	100601 308032	Anti-Histoplasma		1 ml
Meridian Bioscience	100401 308034	Coccidioides "F" Antigen		1 ml
Meridian Bioscience	100801 308035	Anti-Coccidioides "F"		1 ml
Meridian Bioscience	100501 308036	Aspergillus Antigen		1 ml
Meridian Bioscience	100901 308037	Anti-Aspergillus		1 ml
Meridian Bioscience	103001 308038	Coccidioides "TP" Antigen		1 ml
Meridian Bioscience	103101 308039	Anti-Coccidioides "TP"		1 ml
Meridian Bioscience	101009 308040	Immunodiffusion Plates	4 series/plate	9 plates
Meridian Bioscience	101012	Immunodiffusion Plates	1 series/plate	12 plates

Reagents

Legionella

Meridian Bioscience	4873151	FITC-labelled Anti-Human Immunoglobulin	Lyophilized (for 2,5 ml)	1 vial
------------------------	---------	--	--------------------------	--------

Reagents

Meridian Bioscience	XX-715 309001	GULLSORB	IgG inactivation reagent, ready to use	5 ml
Meridian Bioscience	4870141 309002	Mounting medium		3 ml
Meridian Bioscience	4874001 309003	Specimen diluent	Lyophilized (20 ml diluted reagent), recommended for EBV-EA and L. pneumophila IFA	1 vial
Meridian Bioscience	4873691 309004	PBS buffer	Concentrate (30x for 1L working solution)	33 ml